

Tube attachment devices from Hollister provide an easy-to-use alternative to tape for securing tubes, drains, and catheters. Each device features a skin-friendly FlexWear skin barrier to help protect patient skin and an adjustable clamping mechanism to allow for easy tube positioning and better access for site care.



Stock No	Description	Quantity/Box
9781	Horizontal Tube Attachment Device (HTAD)* <ul style="list-style-type: none"> • Secures drains, tubes, or catheters (5-40 Fr) parallel to patient's skin • Secures tubes at exit site or distal to exit site • Sterile, latex-free 	5



Stock No	Description	Quantity/Box
9782	Vertical Tube Attachment Device (VTAD)* <ul style="list-style-type: none"> • Secures drains, tubes, or catheters (5-40 Fr) vertical to patient's skin • Secures tubes at exit site • Sterile, latex-free 	5



Stock No	Description	Quantity/Box
9780	Drain Tube Attachment Device (DTAD) <ul style="list-style-type: none"> • Secures a variety of drains or tubes (10-24 Fr) vertical to patient's skin • Sterile, Latex-free 	5



Stock No	Description	Quantity/Box
9786	Feeding Tube Attachment Device (FTAD) <ul style="list-style-type: none"> • Secures nasogastric/nasointestinal tubes (5-18 Fr) • Butterfly-shaped skin barrier conforms to patient's nose • Unique swivel clamp eliminates pressure on the nares • Latex-free 	20



Stock No	Description	Quantity/Box
9785	Suction Tube Attachment (STAD) <ul style="list-style-type: none"> • Secures nasogastric suction tubes (12-18 Fr) • Butterfly-shaped skin barrier conforms to patient's nose • Unique swivel clamp eliminates pressure on the nares • Non sterile 	20

*Visit <http://www.hollister.com/us/tube/learning/theatre.asp> to view video instructions